

VSV forward primer mix	100 μl	500 μl
Rnase-Free Water	70 μ l	350 μ l
7230F-1 [100 μ M] VSNJV	10 μ l	50 μ l
7230F-2 [100mM] VSNJV	10 μ l	50 μ l
7230F [100mM] VSIV	10 μ l	50 μ l
Total volume	100 μ l	500 μ l

VSV reverse primer mix	100 μl	500 μl
Rnase-Free Water	50 μ l	250 μ l
REV2 [100 μ M] VSNJV	40 μ l	200 μ l
7456R [100mM] VSIV	10 μ l	50 μ l
Total volume	100 μ l	500 μ l

VSV combo probe mix	100 μl	500 μl
Rnase-Free Water	85 μ l	425 μ l
Short [100 μ M] VSNJV	5 μ l	25 μ l
M1 [100mM] VSNJV	2.5 μ l	12.5 μ l
M2 [100mM] VSNJV	2.5 μ l	12.5 μ l
IN 22 [100mM] VSIV	5 μ l	25 μ l
Total volume	100 μ l	500 μ l